

Pacific Southwest Lab Equipment, Inc.



Special Stains Kits

Histology Staining Kits	Use	Contents	Volume	# of tests	Cat. #	Price
<i>A.F.B.</i>	<i>T. bacillus</i>	<i>1% Acid Alcohol Carbol Fuchsin Solution (Store at room temperature) Light Green Counterstain</i>	<i>1 pint 100ml 100ml</i>	100	PSL0011	\$166.00
<i>Alcian Blue, pH 1.0</i>	<i>Mucosubstances</i>	<i>Alcian Blue Stain pH 1.0 Nuclear Fast red stain 3%Acetic Acid</i>	<i>100 ml 100ml 1 pint</i>	100	PSL0022	\$272.00
<i>Alcian Blue, pH 2.5</i>	<i>Mucosubstances</i>	<i>Alcian Blue Stain ph 2.5 Nuclear Fast Red Stain 3% Acetic Acid</i>	<i>100ml 100ml 1pint</i>	100	PSL0033	\$261.00
<i>C.E.M.</i>	<i>Eosinophil</i>	<i>Astra Blue Solution Modified Mayer's Hematoxylin Red Soutlion</i>	<i>100ml 100ml 100ml</i>	100	PSL0044	\$199.00
<i>Colloidal Iron</i>	<i>Mucosubstances</i>	<i>12% Glacial Acetic Acid 3% Hydrochloric Acid Solution 3% Potassium Ferrocyanide Solution Colloidal Iron Stock Solution Van Gieson's Stain</i>	<i>1 pint 100ml 100ml 100ml 100ml</i>	100	PSL0077	\$195.00
<i>Congo Red</i>	<i>Amyloids</i>	<i>Congo Red Amyloid Stain Modified Mayer's Hematoxylin 1% Sodium Hydroxide</i>	<i>100ml 100ml 1 Pint</i>	100	PSL0088	\$192.00
<i>Dieterle's</i>	<i>Spirochete</i>	<i>Wiley's Solution Solution A Hydroquinone 10% Zinc Formalin Alcoholic Gummastic Solution B</i>	<i>100ml 100ml 1.5gm. 100ml 125ml 5ml</i>	100	PSL0100	\$233.00
<i>Elastic</i>	<i>Elastic fibers</i>	<i>5% Alcoholic Hematoxylin Solution Iodine Solution Van Gieson's Stain 10% Ferric Chloride Solution 5% Sodium Thilosulfate Solution</i>	<i>100ml 100ml 100ml 100ml 100ml</i>	100	PSL0122	\$193.00
<i>FITE</i>	<i>Bacteria</i>	<i>Xylene-Peanut Oil 1%ACID Alcohol Carbol Fuchsin Solution (Store at Room Temperature) Methylene Blue Stain</i>	<i>100 ml 1 pint 100 ml 100 ml.</i>	100	PSL0133	\$184.00
<i>Fontana-Masson</i>	<i>Argentaffin cells</i>	<i>0.1 % Gold Chloride Nuclear Fast Red Stain 5% Sodium Thiosulfate Fontana Silver Solution</i>	<i>100 ml 100ml 100ml 5 Vials</i>	100	PSL0144	\$218.00

G.M.S.	Fungi	5% Chromic Acid Solution 0.3% Gold Chloride Light Green Counterstain Methenamine/Silver Compound 1% Sodium Bisulfite 2% Sodium Thiosulfate Sodium Borate Concentrate Slide Jar for Microwave Procedure	100ml 100ml 100ml 5 Vails 100ml 100ml 30ml 4 ea.	100	PSL0155	\$216.00
Giemsa	Bone marrow	May-Grunwald Stock Solution Phosphate Buffer 6.8pH Solution Giemsa Stock Solution	1 pint 1 pint 1 pint	100	PSL0166	\$218.00
Gram McDonald's	Bacteria	Gentian Violet (Primary Stain) Decolorizer Solution Tratrazine Solution Iodine Solution (Mordant) Carbol Fuchsin Solution	100 ml 100ml 100ml 100ml 100ml	100	PSL0188	\$216.00
HELICO* STAT	H. pylori	Periodic Acid Solution Alcian Yellow Solution Sodium Hydroxide Solution Sodium Metabisulfite Solution Toluidine Blue Stock Toludine Blue Working Solution: Distilled Water Toludine Blue Stock	1 Pint 100ml 10ml 100ml. 10ml — 50.0ml 10 drops 3 drops	100	PSL0199	\$311.00
Hepatitis B	Hepatitis	5% Potassium Permanganate Solution 2% Oxalic Acid Solution Differentiation Solution 3% Sulfuric Acid Solution Orcein Stain Solution	30ml 100ml 100ml 30ml 100ml	100	PSL0200	\$272.00
Iron Stain	Ferric Iron	3%Potassium Ferro Cyanide 3% Hydrochloric Acid Nuclear Fast Red Stain	1 pint 1 Pint 100ml	100	PSL0211	\$184.00
Luxol Fast Blue	Myelin sheath	Luxol Fast Blue Stain 70% Reagent Alcohol Cresyl Echt Violet Stain 0.05% Lithium Carbonate	100 ml 1 pint 100ml 1 pint	100	PSL0222	\$276.00
Methyl Green Pyronin	RNA, DNA	Methyl Green Pyronin Procedure	100 ml	100	PSL0233	\$112.00
Mucicarmine	Mucin	Mucicarmine Stock Solution Modified Mayer's Hematoxylin Solution Mentail Yellow Stain Solution	100 ml 100 ml 100ml	100	PSL0244	\$196.00
Oil Red O	Neutral fats	Propylene Glycol Solution Modified Mayer's Hematoxylin Solution Oil Red O Stain Solution (Heat to 60 degrees c prior to use)	1 Pint 100ml 100ml	100	PSL0255	\$188.00
P.A.S.	Glycogen	0.5% Periodic Acid Solution Light Green Counterstain Schief's Solution Modified Mayer's Hematoxylin Solution	100 ml 100ml 1 pint 100ml	100	PSL0277	\$188.00
Von Kossa	Calcium	Nuclear Fast Red Counterstain 5% Sodium Thiosulfate Solution 5% Silver Nitrate Staining Solution	100ml 100ml 100ml	100	PSL0288	\$216.00
P.T.A.H.	Muscle, fibrin	Ferric Ammonium Sulfate Solution P.T.A.H. Stain 10% Zinc Chloride Solution	100 ml 100ml 1 pint	100	PSL0299	\$237.00

<i>Pneumocystis</i>	<i>P. carinii</i>	<i>Cresyl Echt Violet</i> <i>Slide Jar</i> <i>Naphthol Yellow S</i>	100 ml 1 Slide 100ml	100	PSL0300	\$172.00
<i>Reticulum Chandler's</i>	<i>Reticulum fibers</i>	1% Potassium Permanganate Solution 3% Ferric Ammonium Sulfate Solution 0.1% Gold Chloride Solution 20% Formalin 3% Potassium Metabisulfite Solution Nuclear Fast Red Stain 3% Sodium Hydroxide Solution 5% Sodium Thiosulfate Solution	100 ml 100 ml 100ml 100ml 100ml 100ml 100 ml 100 ml 100 ml 5 Vials	100	PSL0311	\$234.00
<i>Sudan Black B</i>	<i>Lipids</i>	Fixative Solution 70% Reagent Alcohol Modified Sudan Black B Solution (Filter as necessary) Modified Mayer's Hematoxylin Solution	100 ml 1 pint 100ml 100ml	100	PSL0322	\$272.00
<i>Copper Stain</i>	<i>Copper</i>	Rhodanine Stock Solution Modified Mayer's Hematoxylin Solution Acetate Buffer Solution	25 ml 100 ml 2 pints	100	PSL0333	\$233.00
<i>Trichome Masson's</i>	<i>Collagen</i>	Bouin's Fluid Aniline Blue Stain Weigert's Hematoxylin "A" Biebrich Scarlet - Acid Fuchsin 1% Acetic Acid Weigert's Hematoxylin "B" Phosphomolybdic/Phosphotungstic Acid	100 ml 100 ml 100 ml 100 ml 100 ml 100 ml 100 ml	100	PSL0344	\$222.00
<i>Steiner</i>	<i>Bacteria</i>	Oxidizer Solution Oxidizer Solution Hydroquinone 10% Zinc Formalin 1% Silver Nitrate	100 ml 100 ml 1.5 gm 100 ml 100 ml 5 ml	100	PSL0355	\$221.00
<i>Warthin-Starry</i>	<i>Bacteria</i>	Spirochete solution (do not refrigerate) 4% Gelatin Solution (refrigerate) 0.1% Hydroquinone, solution (refrigerate) 0.5% Silver Nitrate, solution (refrigerate) 2% Silver Nitrate, solution (refrigerate)	100 ml 75 ml 45 ml 100ml 30 ml	100	PSL0366	\$234.00

Slide Staining Jars call for availability